


# Flexibility and performance of modular “ioniTOF” Time-of-Flight platform used in PTR-TOF and APi-TOF systems

C. Lindinger<sup>1</sup>, L. Märk<sup>1</sup>, J. Herbig<sup>1</sup>, A. Jordan<sup>1</sup>, S. Feil<sup>1</sup>, P. Mutschlechner<sup>1</sup>, G. Hanel<sup>1</sup>, E. Hartungen<sup>1</sup>, S. Jürschik<sup>1</sup>, P. Sulzer<sup>1</sup>.

<sup>1</sup> IONICON Analytik GmbH., Eduard-Bodem-Gasse 3, 6020 Innsbruck, AUSTRIA.





## Proton-Transfer-Reaction - Mass Spectrometry



### ► PTR: Proton Transfer Reaction

- “Soft” ionization of VOCs
  - low fragmentation
  - simple quantification
- Real-time < 100 ms
- Linearity over 6 orders of magnitude
- High sensitivity, LOD < 1 pptv




### ► MS: Time-of-Flight Mass Spectrometry

- High mass resolution
- Separation of isobars
- Identification of chemical composition

## Modular Design – Fit for Flexibility

Time-Of-Flight (TOF) mass spectrometers are state of the art in fields of application where whole mass spectra are needed with high time and mass resolution [1]. Here we present a novel approach of an in-house developed modular TOF platform (ioniTOF) which can be coupled with various sources and inlet systems, e.g. with an established Proton-Transfer-Reaction (PTR) ion source or an Atmospheric Pressure interface (APi).



- Compact orthogonal TOF analyzer with a length of 25 cm
- Ion mirror doubles the flight path → high mass resolution despite small dimensions
- Module of about 20 cm length → Extension of the flight path and increase of the mass resolution

## Ion Funnel and Ion Guide

- Ion-funnel technology: ION BOOSTER
- Series of ring electrodes with decreasing orifices and applied RF field
- Combination with a Hexapole as ION GUIDE is possible
- Modular design → easy upgrading of existing PTR-MS instruments
- Focusing of ions and improvement of ion transmission
- In-house built TOF platform



## Performance Comparison




Figure 1: The sensitivities (cps/ppbv) are obtained by analysis of compounds from a certified gas standard (TO-14A).




Figure 2: Limit-of-Detection (LoD) is calculated using the 3σ (standard deviation) method for integration times between 1 and 60 s.





Figure 3: Left: Mass resolution for protonated Trichlorobenzene (m/z 181). The resolutions are calculated by the FWHM method. Right: Nosespace data of a coffee study demonstrating the high mass resolution of the PTR-TOF 6000 X2. In this example isobaric compounds in the coffee flavor can be clearly separated.





	Funnel	Ion-Guide	TOF
PTR-TOF 1000			compact
PTR-TOF 1000 <i>ultra</i>	✓		compact
PTR-TOF 4000		✓	high resolution
PTR-TOF 6000 X2	✓	✓	very high resolution

## Straight forward ion chemistry


How does the RF field of the ion funnel influence the ion chemistry during the ionization process?

We analyzed the fragmentation pattern of compounds whose product ion branching ratios strongly depend on the reduced electric field in the drift tube (E/N). Two such compounds, octanal and limonene, were measured with and without ion funnel (see Figure below). By adjusting the RF amplitude and the DC voltages of the ion funnel in the PTR-TOF 1000 *ultra* and the PTR-TOF 6000 X2 we could reproduce the branching ratios with an error well below 10%. Our modular ion funnel clearly improves the instrumental sensitivity while not considerably changing the ion chemistry. Branching ratios and calculated concentrations remain accurate.



## APi-TOF

We have developed a novel Atmospheric Pressure interface (APi; see scheme below) for the analysis of positive as well as negative ions in the environment. The system, which has been designed to be coupled with the ioniTOF platform, consists of an atmospheric pressure interface with a critical orifice for contact free sample introduction and a hexapole ion guide for high ion transmission efficiency over a broad m/z range (in contrast to common quadrupole ion guides). By combining two virtually identical versions of this setup with opposite polarity into one rugged and compact instrument, both ion polarities can be detected simultaneously in real-time.



### References

[1] A.M. Ellis, C.A. Mayhew, In: Proton Transfer Reaction Mass Spectrometry: Principles and Applications, Chichester, UK: John Wiley & Sons, Ltd (2014). Basisprogramm and an ASAP project.

### Acknowledgement

This project was supported by the Land Tirol via a FuEul project and by the Austrian Research Promotion Agency (FFG) via a