

Characterization of a Novel Laminar-Flow Oxidation Reactor for Simulating Atmospheric Multiple-Day Oxidation


Tobias Reinecke, Markus Leiminger, Andreas Klinger, Martin Graus and Markus Müller

IONICON Analytik GmbH., Eduard-Bodem-Gasse 3, 6020 Innsbruck, AUSTRIA




IONICON


Atmospheric Multiple-Day Oxidation



FUSION PTR-TOF 10




CHARON Particle Inlet



IONICON Laminar-Flow Oxidation Reactor – ILOx




The automatable IONICON Laminar-flow Oxidation Reactor (ILOx) consists of a 110 cm long Quartz glass tube with a total internal volume of 8 l. Typical residence times range from 2 to 16 min. All materials are optimized for purest experiments possible (Quartz glass, Teflon and passivated surfaces only). The flow reactor is surrounded with a powerful 300 W UVA LED lamp (365 nm) for direct irradiation over 80 cm of the reactor. To maintain stable temperatures inside the ILOx reactor active cooling via a set of fans is fully sufficient. To achieve a laminar flow condition inside the reactor, the carrier gas is injected by a circular motion, suppressing any formation of injection jets. Various oxidants can be supplied via multiple customizable inlet ports, e.g. ozone can be generated externally from O_2 by UVC irradiation.



Flow Conditions

The geometries for optimal sample injection and extraction were optimized via 3D computational fluid dynamics (CFD) simulations using OpenFOAM (openfoam.org).



> Inlet

Our 3D CFD simulations resulted in a best initial mixing and following optimal laminar flow conditions via a dual radial injection of the buffer gas and the sample. A benefit of this radial injection is the prevention of a formation of a jet at the center of the reactor.


> Reactor

We have chosen a length-to-diameter ratio of > 10 for reaction conditions that resemble a traditional aerosol-flow-tube reactor. With the optimized inlet geometry optimal laminar flow conditions can be reached at various residence times.


> Outlet

A conical outlet is followed by a core sampling of VOCs to minimize the loss of low volatiles while simultaneously transmitting the particles at highest efficiencies.

Laminar-Flow Characterization



Particle Transmission Efficiency




> Measurement of Ammonium Sulfate (AS) particles before and after ILOx.


> Excellent transmission of AS particles almost independent of particle size.

Rapid Oxidation of Xylene

- > 8 ppbV of o-xylene were injected to study the rapid oxidation potential of the ILOx reactor.
- > IPN was added as a OH precursor. With active UVA irradiation IPN is photolyzed to produce OH radicals and acetone.
- > This allows for a rapid oxidation of xylene resembling an atmospheric oxidation of approximately two days.
- > MCM 3.3.1 simulations of this experiment result in average OH concentrations of 3×10^8 cm⁻³, which equals to an OH exposure OH_{exp} of 10^{11} cm⁻³s.



Modelling of xylene decay and radical formation (MCM 3.3.1)



References

- [1] Tobias Reinecke, Markus Leiminger, Alfons Jordan, Armin Wischler, and Markus Müller, 2023, Ultrahigh Sensitivity PTR-MS Instrument with a Well-Defined Ion Chemistry, Anal. Chem., 95, 32, 11879–11884, <https://doi.org/10.1021/acs.analchem.3c02669>
- [2] Tobias Reinecke, Markus Leiminger, Andreas Klinger, and Markus Müller, Direct detection of condensed particulate polycyclic aromatic hydrocarbons on a molecular composition level at low pg m⁻³ mass concentrations via proton-transfer-reaction mass-spectrometry, Aerosol Research, 2, 225–233, 2024.