

# FUSION PTR-TOF 10k – Next generation PTR-TOF instrument



IONICON Analytik Ges.m.b.H.

Eduard-Bodem-Gasse 3, 6020 Innsbruck, AUSTRIA

<https://www.ionicon.com/products/details/fusion-ptr-tof-10k>


## Introduction

The possibility of on-line measurement and detection of compounds with high sensitivity has been highlighted as one of the main advantages of Proton-Transfer-Reaction – Mass Spectrometry (PTR-MS) [1] since its introduction in the 1990s.

The quest for the highest sensitivity has recently led to a multitude of reaction chamber designs, all of which have in common that they employ RF fields to guide and focus the ions. However, this comes with a major drawback, namely that the well-defined ion chemistry in PTR-MS, which is one of the most important features of this technology, is sacrificed [1].

Here we present the FUSION PTR-TOF 10k, which is part of IONICON's Next-Gen PTR-TOF instrument series and combines the novel TRION source and FUSION drift cell. The novel TRION source plus the FUSION reaction chamber and a high resolution TOF analyzer are the new benchmark for nearly instant reagent ion switching, ppqv limit of detection and reliable E/N-values.

## FUSION PTR-TOF 10k



### FUSION reaction chamber

- ⊗ mounted within an enveloping recipient
- ⊗ directly injected sample air flows through a series of gas-tight ring electrodes
- ⊗ pumped through the gaps between the adjacent ion funnel electrodes
- ⊗ leaks are virtually impossible due to enveloping recipient

### TRION source


- ⊗ 3 reagent ion sources, positioned around a virtually contact free sample inlet
- ⊗ gas flow in direction of drift tube
- ⊗ continuously producing up to 3 reagent ions at the same time, e.g.  $\text{H}_3\text{O}^+$  [1],  $\text{NO}^+$  [2],  $\text{NH}_4^+$  [3],  $\text{O}_2^+$  [2]
- ⊗ switching times up to 100-200ms

## Key Advantages

### A unique combination of:

- ⊗ high sensitivities up to 30,000 cps/ppbv
- ⊗ LoDs down to <100 ppqv
- ⊗ selective reagent ion (SRI) switching in about 200ms, high stability of different reagent ion sources
- ⊗ FUSION PTR technology is best combined with ionTOF 10k with a mass resolution of typically 10,000 – 15,000  $m/\Delta m$  (FWHM)
- ⊗ patented TRU-E/N ion chemistry quality

## Overview



### A Next-Gen PTR-TOF setup


- ⊗ ultra-clean FUSION reaction chamber
- ⊗ TRION source for instant SRI switching, stable ionization conditions due to continuous reagent ion production
- ⊗ advanced inlet design for milliseconds response time, gas flow in direction of drift tube
- ⊗ TRU-E/N ion chemistry quality

## Results

The performance of the novel FUSION PTR-TOF 10k instrument was evaluated with a gas feed containing 1-2 ppbv of various compounds.


### Fast selective reagent ion (SRI) switching

With the new setup, reagent ions can be switched with times as short as 200ms, as shown below for switching the reagent ion  $\text{H}_3\text{O}^+$  to  $\text{O}_2^+$ .



### Mass resolution and high sensitivities

At the same time, high mass resolution and high sensitivities are achieved time, as shown below for  $\text{C}_6\text{H}_3\text{Cl}_3\text{H}^+$   
Mass resolution: >12,000 (FWHM)  
Sensitivity: > 10,000 cps/ppbv



### TRU-E/N

We were able to reproduce the fragmentation pattern of pinonic acid upon PTR ionization. The branching ratios are very sensitive to the applied E/N. Utilizing our patented TRU-E/N method (US10074531, EP3309817), we were able to reproduce intensity ratios for all relevant pinonic acid product ions corresponding to E/N values from regular operating conditions, i.e. without the use of RF voltages, of 120 Td down to 60 Td using a FUSION PTR-TOF.

For example, using the relative ion yield of the thirteen different fragments formed from pinonic acid during ionization by PTR, we obtain a Pearson correlation coefficient of 0.97 in the case of 60 Td.

### References

- [1] A.M. Ellis, C.A. Mayhew, In: Proton Transfer Reaction Mass Spectrometry: Principles and Applications, Chichester, UK: John Wiley & Sons, Ltd (2014). ISBN: 978-1-405-17668-2
- [2] A. Jordan et al., An online ultra-high sensitivity Proton-transfer-reaction mass spectrometer combined with switchable reagent ion capability (PTR + SRI – MS). Int. J. Mass Spectrom. 286 (2009) 32-38. DOI: 10.1016/j.ijms.2009.06.006
- [3] M. Müller et al., A novel method for producing  $\text{NH}_4^+$  reagent ions in the hollow cathode glow discharge ion source of PTR-MS instruments. Int. J. Mass Spectrom. 447 (2020) 116254. DOI: 10.1016/j.ijms.2019.116254