


# Advantages of Proton Transfer Reaction – Mass Spectrometry (PTR-MS) in the Analysis of Potentially Dangerous Substances

P. Sulzer<sup>1</sup>, A. Jordan<sup>1</sup>, G. Hanel<sup>1</sup>, E. Hartungen<sup>1</sup>, H. Seehauser<sup>1</sup>, L. Märk<sup>1</sup>, S. Haidacher<sup>1</sup>, R. Schottkowsky<sup>1</sup> and T. D. Märk<sup>1,2</sup>

<sup>1</sup> IONICON Analytik Gesellschaft m.b.H., Technikerstr. 21a, A-6020 Innsbruck, Austria

<sup>2</sup> Institut für Ionenphysik und Angewandte Physik, Leopold-Franzens Universität, Technikerstr. 25, A-6020 Innsbruck, Austria




## Abstract

Proton Transfer Reaction – Mass Spectrometry (PTR-MS) is an innovative technology invented and developed by scientists of the "Institut für Ionenphysik" at the Leopold-Franzens University in Innsbruck and made commercially available by the spin-off company IONICON Analytik GmbH. Details about the technique have been described elsewhere (e.g. [1], [2] and very recently [3]). In short, in a hollow cathode ion source hydronium ions (in the novel "SRI" mode also NO<sup>+</sup> or O<sub>2</sub><sup>+</sup>) are produced at very high purity levels (up to 99.5%) and afterwards injected into a drift tube where the actual ionization process takes place via proton transfer from the hydronium (or charge transfer from NO<sup>+</sup> or O<sub>2</sub><sup>+</sup>) to the trace constituents. Finally either a quadrupole or a time-of-flight mass spectrometer analyzes the product ions according to their masses and yields. With this setup trace gas compounds can be measured in a concentration range from several ppmv down to the ppqv (parts-per-quadrillion) region [4] with a typical response time well below 100ms and, in case a TOF mass analyzer is used, a mass resolution better than 6.000 m/Δm. One can easily see that these performance data would make the PTR-MS technology an ideal candidate for the detection of illicit substances such as explosives, CWAs, drugs. Here we present first results of measurements on explosives and a chemical warfare (CWA) analogue utilizing a high mass resolution and high sensitivity PTR-TOF-MS instrument as well as a quadrupole based PTR-MS instrument. In addition we present data obtained with our novel "inlet booster system", which improves the detection limit by a factor of 500 already in its preliminary version.




### Inlet booster system

The sample air is drawn through a metal mesh at a high flow rate for about 10s. Afterwards the mesh is resistively heated and the evaporating compounds are analyzed with a common PTR-MS instrument. The whole process takes less than one minute and can therefore still be called "online".




### Direct inlet

Direct inlet measurement on a high mass resolution time-of-flight (PTR-TOF-MS) instrument. TNT on protonated mass 228.025m/z can be clearly separated from an unknown lab-air compound.



### Inlet booster



Performance check of the inlet booster system (PTR-Quad-MS). The left diagram shows a measurement on TNT, at first via the direct inlet (~200cps), then a check if the booster is clean and finally a 55s boosting procedure that leads to an increase in signal intensity by a factor of ~500. The right diagram demonstrates that the principle also works for other explosives (Note: the booster was not cleaned between the three samples).



Mass spectrum of TNT utilizing a Compact PTR-Quad-MS (reduced in size and weight but also in sensitivity). The protonated parent ion on mass 228m/z can still be detected. On mass 183m/z we observe protonated DNT, which is a known impurity of commercial TNT and has a significantly higher vapor pressure.

Tests on a few mg of TNT in a glass vial utilizing a PTR-TOF-MS. Top left: direct inlet starting at about cycle 50; top right: utilizing the inlet booster system on open vial; bottom left: utilizing inlet booster system on closed vial; bottom right: utilizing inlet booster system on human finger that was in contact with TNT (first ~100 cycles: check if booster is clean).

Part of a mass spectrum obtained with a PTR-Quad-MS of a vial that was flushed with C<sub>4</sub>H<sub>8</sub>ClS (compare: mustard gas C<sub>4</sub>H<sub>8</sub>Cl<sub>2</sub>S). Note the extremely high signal intensity and the nicely matching isotope ratio between mass 125 and 127m/z that makes an unambiguous identification possible.



### References

- [1] W. Lindinger, A. Hansel, A. Jordan, Int. J. of Mass Spectrom. and Ion Processes; 173/3 (1998) 191-241.
- [2] J. de Gouw, C. Warneke, T. Karl, G. Eerdeken, C. van der Veen, R. Fall, Mass Spectrometry Reviews, 26 (2007), 223-257.
- [3] R. S. Blake, P. S. Monks, A. M. Ellis: Proton-Transfer Reaction Mass Spectrometry, Chem. Rev., 109 (3) (2009), 861-896.
- [4] A. Jordan, S. Haidacher, G. Hanel, E. Hartungen, L. Märk, H. Seehauser, R. Schottkowsky, P. Sulzer and T.D. Märk, Int. J. Mass Spectrom., 286 (2009) 32-38.